

ATTENTION ALL CLINICAL STAFF (MD, RN, APP)

DO YOU HAVE QUESTIONS FOR RESPIRATORY THERAPISTS REGARDING THE USE OF YOUR VENTILATORS OR EQUIPMENT?

STEP 1- CALL THE CLINICAL COMMAND CENTER

- Call the Clinical Command Center at **1-800-TO-SINAI (867-4624)**

STEP 2 – ASK FOR TELE-RESPIRATORY SUPPORT HOTLINE

- Ask for the **Tele-Respiratory Support Hotline**
- You will be connected to the MSHS On-Call Respiratory Therapist Team

STEP 3 – LOCATE YOUR InTouch HEALTH TELEMEDICINE DEVICE

- Note: There are two devices on your unit that you can use:
 - InTouch Light 4 device
 - Apple iPad on a stand


InTouch Lite 4 Device



STEP 4 – LOCATE THE NAME OF YOUR DEVICE

- Locate the name of your InTouch device (label on the back of the device).
- Here are the naming conventions:
 - Lite-4: Facility Lite-4 Location #
 - *Ex. MSB Lite-4 ED 1*
 - iPad: Facility Location Patient #
 - *Ex. MSH KP 6E Med Surg Patient 1*

Apple iPad



STEP 4 – BRING InTouch DEVICE TO PATIENT

- Bring your InTouch device to where the patient is, in preparation for the Tele-Respiratory support call.

STEP 5 – TELL THE HOTLINE THE NAME OF YOUR DEVICE

- Tell the hotline support person the name of your device.
- The respiratory therapist will beam into your device to assist you.